

Unsurpassed Performance in Medical Radiography

The Choice for Versatility, Reliability and Value



AM3 Integrated Table System

AMERICOMP

Affordable Choices for General Radiography

Since 1987, the quality standards of Americomp have been unsurpassed. Americomp systems remain the affordable choice in today's marketplace, offering the same long-term dependable performance that has become synonymous with the Americomp name. Our manufacturer's parts warranty secures your investment today and guarantees your investment value tomorrow.



AM3 Integrated Table System Features:

- Tubestand integration to table saves valuable space and does not require wall or ceiling modifications
- Unique off-table option provides a full range of off-table procedures including standing knee and ankle exams (shown above)
- Versatile tube positioning offers a wide range of procedures
- Flat, extra-wide, table top allows easy patient positioning
- Heavy-duty construction with 300 lb. patient rating
- Wallstand provides effortless positioning of stationary grid cabinet or optional reciprocating bucky



AM4 Stationary Table System Features:

- Floor-to-wall or ceiling tubestand offers versatile tube positioning including off-table standing knee and ankle exams
- Flat, extra-wide, table top allows easy patient positioning
- Heavy-duty construction with 300 lb. patient rating
- Wallstand provides effortless positioning of stationary grid cabinet or optional reciprocating bucky



AM5 Compact Table System Features:

- Perfect for "no space" situations
- Fits into rooms measuring as small as 7 ft. x 9 ft.
- Floor-to-wall or ceiling tubestand offers versatile tube positioning
- Flat table top allows easy patient positioning
- Heavy-duty construction with 300 lb. patient rating
- Wall mounted cassette holder saves space and does not require handling of cassette tray

High Performance Generators

Complement any Americomp radiographic system with a high-quality, state-of-the-art generator. Contact your Americomp dealer for details.

Digital-ready Systems

Americomp radiographic equipment is completely compatible with most digital imaging systems. Digital imaging technology may be added at any time.



© 2020 SMIC. All rights reserved.
1721 Stewart Street,
Santa Monica, CA 90404
Toll Free: (800) 859-6162
Email: info@spectrumxray.com
Website: www.spectrummedimaging.com